



HONG KONG
ASIAWORLD-EXPO
亞洲國際博覽館

3RD TO 6TH
DECEMBER
2025






Innovation in Reprocessing Robotic Surgical Instruments

Name: Nick Dassatti
Affiliation: Intuitive



Evolution of Robotic Surgery and Reprocessing



Manual cleaning

Automated cleaning for
robotic instruments

Pre-cleaning
devices



Robotic Surgical Innovation Drives Reprocessing Innovation

Development of **new solutions** to **improve** robotic reprocessing and sterile processing workflows.

- Automated ultrasonics with flushing and articulation
- Pre-cleaning pumps
- Drying cabinets
- Protection and storage




Reprocessing Innovations and Patient Safety

- Validation of state-of-the-art reprocessing solutions **demonstrates high level of patient safety.**
- 10 years of PQ data on robotic instruments **resulted in an average of 34µg total protein per device** using validated washer-disinfectors.
- New manual pre-cleaning technologies **show improvement of reprocessing workflow** with reduced cleaning steps.



Example of PQ testing




Cleaning of robotic instruments:
Can we reduce the work load in
the CSSD and improve patient
safety?

Klaus Roth; Henri Hubert; Sylvia Pfeiffer; Silke Buchmann; Ludger Schrieder

SMP GmbH Prüfen Validieren Forschen
Hechingen Strasse 262
72072 Tübingen; Germany




Historical Approach to Development



Early design iterations



Instrument manufacturer design freeze



Reprocessing solutions

Limitations

- Requires years of development and validation.
- Limited by timing of final device design.
- Dedicated robotic racks and washer-disinfectors.
- Specialized cleaning cycles.



Development of a New Approach to Automated Cleaning

Can we improve on existing solutions and the current development process?

Goals

- Maintain or improve on the high level of patient safety previously established.
- Improve workflow and equipment compatibility for users.
- Faster time-to-market of reprocessing solutions to support site case load.

How do we get there?

Key Questions

- Can we achieve the high bar set by existing reprocessing solutions?
- What equipment do global users have?
- How do we accommodate for the diversity of washer-disinfectors early into the design process?



Development of a New Approach to Automated Cleaning

How do we get there?

Key Questions

- Can we achieve the high bar set by existing reprocessing solutions?
- What equipment do global users have?
- How do we accommodate for the diversity of washer-disinfectors early into the design process?



Approach


- Comprehensive feasibility testing and formal validation to international standards.
- In depth review of equipment common among users and sites with robotics.
- Compliance to ISO standards.




Solution – Instrument Reprocessing Trays



Reprocessing Tray Approach for Robotic Instruments



Early design iterations



Automated cleaning parameters per ISO 17664-1 with connection method



Benefits

- More solutions available sooner to the user.
 - Washer-disinfector flexibility based on connection design.
 - Enables routine cycle parameters.
- Less dedicated space to robotic reprocessing.
- Improved efficiency with washer-disinfectors being used for both robotics and non-robotics.
- Multi-purpose with sterilization compatibility.



Validations of Cleaning Efficacy

**Surgical
Simulated Use
and Drying**



**Manual
Pre-Cleaning**




**Automated
Cleaning**



Extraction



Analysis



Residuals analysis
in accordance with
ISO 15883-5: 2021.

Protein:

- $< 6.4 \mu\text{g}/\text{cm}^2$
- $< 200 \mu\text{g}/\text{device}$

Hemoglobin:


- $< 2.2 \mu\text{g}/\text{cm}^2$



Validations of Steam Sterilization Efficacy

Test Method and Results

- Three cycles with fully loaded trays.
- Worst-case biological indicator placement.
- Overkill method for a 6-log reduction in a half-cycle.
- No growth after incubation period.



Innovation Outcomes – PQ Testing

Cleaning Performance Qualification (PQ) Testing

In accordance with:

- HTM 01-01Part D
- ISO 15883-1: 2025 / ISO 15883-5: 2021

Collected data in clinical settings to <100ug levels as seen with historical data for innovative approaches to robotic reprocessing.

UK Quarterly Periodic Testing - Before



UK Quarterly Periodic Testing - After



EU PQ Testing - Extraction

Innovation Outcomes - Implementation

SteriLog – Lucerne, Switzerland

Challenged with limited time to start reprocessing of new robotic instruments.

- Tray approach allowed for fast implementation.
- Simple integration into existing workflow.
- Minimal adjustment to existing wash rack.
- Quick turnaround of successful performance qualification (PQ) testing.



Markus Seichter
Managing Director of SteriLog
Sterile Goods Supply Lucerne AG



Innovation Outcomes – Handling

London, UK

- Trays provide security of instruments throughout cleaning and sterilization.
- Reduction in manual handling of products.
- Alleviates space constraints in a small, busy sterile processing department.



Future Opportunities

- Continued collaboration for collection of PQ data.
- Cross collaboration across companies:
 - Further optimization of cycle parameters.
 - Combination of tray with other pre-cleaning technologies.
- Design to fit future device designs.



Questions?





Disclosures

Intuitive technology should only be used by physicians and staff who have received specific training in the use of Intuitive technology. Training provided by Intuitive is limited to the use of Intuitive technology, instruments, and accessories. It does not replace the necessary medical training and experience required to perform surgical procedures. Nothing in Intuitive's training curricula makes expressed or implied promises or commitments to the persons attending. Some of the training described in Intuitive promotions may be provided by organizations or entities other than Intuitive.

Some products, features or technologies may not be available in all countries. Please contact your local Intuitive representative for product availability in your region. Refer to the product specific User Manual for indications, contraindications, warnings and other product information. In the U.S., for important safety information, indications for use, risks and full cautions and warnings, please refer to the user manual(s) and www.intuitive.com/safety.

Reprocessing

Refer to the manufacturer's reprocessing instructions for full details regarding cleaning and disinfection of the items listed above. It is the responsibility of the owner of the da Vinci surgical system to properly train and supervise its personnel to ensure that the da Vinci instruments, accessories, endoscopes, and cameras are properly cleaned, disinfected, and sterilized as required by the reprocessing instructions. The da Vinci products should not be used in a clinical setting unless the institution has verified that these devices are properly reprocessed in accordance with the da Vinci system reprocessing instructions.

Third-party products

Refer to your da Vinci surgical system, instruments, and accessories User Manuals for information regarding compatibility of third-party products. Not all third-party products selected for use by the surgeon have been validated or are endorsed by Intuitive for use with the da Vinci system. Some third-party products selected for use by the surgeon may not be available in all countries.

Intuitive's Privacy Notice is available at www.intuitive.com/privacy.

© 2025 Intuitive Surgical Operations, Inc. All rights reserved. Product and brand names/logos, including Intuitive, Da Vinci, and Ion, are trademarks or registered trademarks of Intuitive Surgical or their respective owner.

