



**HONG KONG**  
ASIAWORLD-EXPO  
亞洲國際博覽館

3<sup>RD</sup> TO 6<sup>TH</sup>  
DECEMBER  
2025



# **Investigating Surgical Instrument Damage: How to ensure your investment is protected**

Name: Matthias Tschoerner

Affiliation: Chemische Fabrik Dr. Weigert GmbH & Co. KG, Hamburg



# Topics




- Types of corrosion
- Instrument damage by corrosion
- Common causes and influencing factors
- Examples from praxis, prevention and correction



# Types of corrosion

- Pitting corrosion
- Wear / fretting corrosion
- Stress corrosion
- Foreign and flash rust / secondary rust
- Crevice corrosion
- Surface corrosion
- Contact corrosion




# Most common types of corrosion

## Type of corrosion, appearance

## Reason and influencing factors


Pitting corrosion      Pinhole like black points with rainbow colored surrounding or rust

Chlorides,  
moisture,  
temperature,  
time,  
Material



Wear / fretting corrosion      Scoring and scratches, most in joints

Construction,  
insufficient care  
with oil,  
Material





# Most common types of corrosion

## Type of corrosion, appearance

Stress corrosion      Stress cracking and fractures, often on joints or proximal end of forceps

Foreign or flash rust / secondary rust      Spot like up to uniform brownish discoloration

Crevice corrosion      Rust within cavities

Surface corrosion      Uniform discoloration

Contact corrosion      Spot like brownish discoloration

## Reason and influencing factors

Material stresses, material, temperature, chlorides



Processing unalloyed steels and / or corroded steels, water rust contaminated

Medium, temperature, time, gap width, material

Medium, temperature, time, material

Medium, different materials in contact



# Pitting corrosion

## Action mode:



- Halides esp. chlorides penetrate the passive layer of stainless steel
- Corrosion starts within minutes

## Sources:


- Isotonic solution (i.e. physiologic saline)
- Blood, saliva, and secretion
- Water for processing

## Influencing factors:


- Temperature
- Time / delayed reprocessing



# Pitting corrosion



# Circulation of instruments





# Test of chloride resistance

Dr. Gerhard Kirmse  
Surface changes - causes, consequences,  
solutions: Presentation at German Society  
Sterile Material Supply, 2015


## Test of stainless-steel instruments

- 2 – 100 mg/l chloride
- Boiling water (thermal disinfection)

ISO 13402:1995 - Surgical and dental hand  
instruments — Determination of resistance  
against autoclaving, corrosion and thermal  
exposure

## Results

- 20 mg/l chloride visible signs  
of corrosion
- 100 mg/L chloride – clear corrosion  
of instruments



# Water for the processing of medical devices

## Washer Disinfector

- Tap water (utility water)
- Softened water (utility water)
- Deionised water (critical water)

ANSI/AAMI ST108:2023 Water for the processing of medical devices

Working group for instrument reprocessing (AKI), [www.a-k-i.org](http://www.a-k-i.org)

## Steam Sterilizer


- Deionised water

ANSI/AAMI ST108:2023

Working group for instrument reprocessing (AKI), [www.a-k-i.org](http://www.a-k-i.org)

EN 285:2015 (A1:2021) Steam sterilizers – Large sterilizers

American  
National  
Standard



# Recommended Water Qualities

Flushing / Washing / Rinsing:

Softened Water (Utility Water)

	AAMI ST108	AKI ( <a href="http://www.a-k-i.org">www.a-k-i.org</a> )
Appearance	Colorless, clear, without sediment	Colourless, clear
pH-value	6.5 – 9.5	5 – 8
Total hardness	< 150 mg CaCO <sub>3</sub> /l	< 53 mg CaCO <sub>3</sub> /l (3 °d)
Chlorides	< 250 mg/l	≤ 100 mg/l
Evaporation residue		≤ 500 mg/l
Conductivity	≤ 500 µS/cm	



# Recommended Water Qualities

## Final Rinse: Deionized water (Critical water)

	AAMI ST108	AKI ( <a href="http://www.a-k-i.org">www.a-k-i.org</a> )
Appearance	Colourless, clear, without sediment	Colourless, clear, without sediment
Conductivity	< 10 $\mu\text{S/cm}$	$\leq 15 \mu\text{S/cm}$
pH-value	5 – 7.5	5 – 7.5
Water hardness	< 1 mg $\text{CaCO}_3/\text{l}$	$\leq 2 \text{ mg CaCO}_3/\text{l}$ (0.02 mmol/l CaO)
Evaporation residue		$\leq 10 \text{ mg/l}$
Chlorides	<b>&lt; 0.1 mg/l</b>	<b><math>\leq 0.5 \text{ mg/l}</math></b>
Silicates ( $\text{SiO}_2$ )	< 1 mg/l	$\leq 1 \text{ mg/l}$
Phosphate ( $\text{P}_2\text{O}_5$ )	< 1 mg/l	$\leq 0.5 \text{ mg/l}$








# Corrosion under the microscope

## Risks from corrosion and instrument damage


- Cleaning efficacy critical
- Sterilization efficacy critical
- Rust transfer and further damage of intact instruments
- Costs for repair and exchange

→ Costs and consequential damage










- **Starting - Heavy contamination**
- **Reason - Delayed reprocessing**
- **Insufficient precleaning / cleaning**

## Orthopedics



# Point of use treatment – Corrosion prevention by instrument foam?



## INSTRUCTIONS\*

\* Check with medical device manufacturers instructions for use and material compatibility


- Spray undiluted directly onto the surgical instruments after surgery
- Ensure complete wetting of instruments





# Use of Foam Spray?

- **Analysis of Stainless Steel and Anodized Aluminium Material Compatibility with Foam Sprays Used for Keeping Used Surgical Instruments Moist**, *Information from the "Working Group Instrument Preparation" (AKI), H. Biering, W. Fuchs, J. Staffeldt, Central Service 2010; 18 (4) : 235–243*




**Fig. 1:** Ring cell placement on test pieces



# Corrosion protection test

## DIN 51360

**Corrosion protection test for the protection of iron or steel-based materials (products, tools, machines) in contact with water-mixed cooling lubricants, drilling emulsions, cleaners, cooling and quenching media.**



# Corrosion protection test





*Conducted by an accredited third-party institute*



## 2 Performance

Process	Corrosion chip filter test according to DIN 51360-2** (1981-07) Double determination
Device	2g dry gray cast iron chips GG 25, Exposure time 2h to 2 ml medium at RT

## 3 Test results

Sample	Result 1. Test	Result 2. Test
1. deionized water (conductivity: 0,4 µS/cm)	 (Net weight chips: 1,98 g)	 (Net weight chips: 2,01 g)
2. Saline solution 0.9% m/V in deionized water	 (Net weight chips: 2,02 g)	 (Net weight chips: 2,02 g)



# Corrosion protection test

*Pure Product Application*



Sample	Result 1. Test	Result 2. Test	Sample	Result 1. Test	Result 2. Test
A	 (Net weight chips: 2,00 g)	 (Net weight chips: 1,98 g)	E	 (Net weight chips: 1,97 g)	 (Net weight chips: 2,01 g)
B ✓	 (Net weight chips: 2,02 g)	 (Net weight chips: 2,00 g)	F	 (Net weight chips: 2,03 g)	 (Net weight chips: 1,98 g)
C	 (Net weight chips: 1,99 g)	 (Net weight chips: 1,97 g)	G ✓	 (Net weight chips: 2,01 g)	 (Net weight chips: 2,01 g)
D	 (Net weight chips: 2,00 g)	 (Net weight chips: 2,02 g)	H	 (Net weight chips: 2,00 g)	 (Net weight chips: 1,99 g)





















# Corrosion protection test

Product + 0.9 % saline solution



Sample	Result 1. Test	Result 2. Test
A	 <p>(Net weight chips: 1,98 g)</p>	 <p>(Net weight chips: 1,99 g)</p>
B	 <p>(Net weight chips: 1,99 g)</p>	 <p>(Net weight chips: 2,02 g)</p>
C	 <p>(Net weight chips: 2,00 g)</p>	 <p>(Net weight chips: 1,98 g)</p>
D	 <p>(Net weight chips: 2,00 g)</p>	 <p>(Net weight chips: 2,02 g)</p>


Sample	Result 1. Test	Result 2. Test
E	 <p>(Net weight chips: 1,98 g)</p>	 <p>(Net weight chips: 1,98 g)</p>
F	 <p>(Net weight chips: 1,99 g)</p>	 <p>(Net weight chips: 2,02 g)</p>
G	 <p>(Net weight chips: 2,02 g)</p>	 <p>(Net weight chips: 2,00 g)</p>
H	 <p>(Net weight chips: 1,99 g)</p>	 <p>(Net weight chips: 2,00 g)</p>



# Prevention - corrosion inhibition by instrument foam?


- ✓ Suitable Pre-wetting in operation rooms
- ✓ Corrosion inhibitor prevents pitting corrosion
- ✓ Prevents corrosion up to 72 hours

Without corrosion inhibition, 0,2 ml NaCl  
7x magnification      40x magnification



With corrosion inhibition 0,2 ml NaCl + 0,2 ml/l,  
7x magnification  
→ pitting corrosion has been avoided 72 h





- **Damage - Corrosion & bone material**
- **Correction - acidic cleaning**
- **Prevention - Revised procedure, fast reprocessing, use of anti-corrosion foam**



# Conclusion

- Corrosion causes hygienic and infection risks
- Corrosion causes costs for repair and exchange
- Common corrosion causes are found in the instrument cycle
- ✓ Contact with isotonic solutions should be avoided
- ✓ Chlorides in the water for reprocessing should be as recommended
- ✓ Instrument foams should be tested for corrosion inhibition
- ✓ Joints (metal-on-metal sliding surfaces) should be lubricated in every reprocessing cycle





Thanks to:  
D. Eisert,  
I. Haacke,  
Dr. B. Wulff

Thank you!



Dr. Matthias Tschoerner  
Chemische Fabrik Dr. Weigert GmbH & Co. KG Hamburg (Germany)

