

25th wfhss CONGRESS



20-23
NOV 2024
SANTIAGO-CHILE

Establishing preconditions for effective duodenoscope reprocessing: an observational cohort study

Koen van der Ploeg^{1,2}, Margreet C. Vos², Nicole S. Erler^{3,4}, Adriana J.C. Bulkman⁵
Bibi C.G.C. Mason-Slingerland², Juliëtte A. Severin^{2*}, Marco J. Bruno^{1*}

1. Dept. of Gastroenterology and Hepatology, Erasmus MC University Medical Center, Rotterdam, The Netherlands.
2. Dept. of Medical Microbiology and Infectious Diseases, Erasmus MC University Medical Center, Rotterdam, The Netherlands.
3. Depart. of Biostatistics, Erasmus MC University Medical Centre, Rotterdam, The Netherlands
4. Depart. of Epidemiology, Erasmus MC University Medical Centre, Rotterdam, The Netherlands
5. Depart. of Quality assurance and regulatory affairs office Medical Technology, Erasmus MC University Medical Centre, Rotterdam, The Netherlands

ERASMUS MC UNIVERSITY MEDICAL CENTRE ROTTERDAM, THE NETHERLANDS




**Adriana J.C. (Diana)
Bulkman**

Expert on sterile medical devices and
endoscope reprocessing (since 2016)

Member of SFERD (since 2018)




THE ISSUE



1. Rauwers et al, High prevalence rate of digestive tract bacteria in duodenoscopes: a nationwide study. Gut. 2018

CURRENT REPROCESSING



GUIDELINES

- Lack of uniformity
- Lack of sufficient evidence
- Essential risk factors remain unidentified
- Lack of information about individual competencies for reprocessing staff




The Digestive Endoscopy Society of Taiwan (DEST)



PROFESSIONAL STANDARD (SFERD)


Kwaliteitshandboek
Reiniging en Desinfectie
Flexibele Endoscopen
Versie 6.0, 2022



Koppel Medische Technologie
SVN
DSMH
VHIG
V&VN

[SFERD Handboek 6.0 NL \(infectiepreventieopleidingen.nl\)](http://infectiepreventieopleidingen.nl)

Professional Standard Handbook
Cleaning and Disinfection
Flexible Endoscopes
Version 4.1, 2017



Koppel Medische Technologie
SVN
DSMH
VHIG
V&VN

[SFERD Handbook 4.1 ENG \(infectiepreventieopleidingen.nl\)](http://infectiepreventieopleidingen.nl)

MANUAL PROFESIONAL
ESTANDAR
ENDOSCOPIOS FLEXIBLES
Limpieza y Desinfección
Versión en español 3.2 2018



Koppel Medische Technologie
SVN
DSMH
VHIG
V&VN

[SFERD Handbook 3.2 ESP \(infectiepreventieopleidingen.nl\)](http://infectiepreventieopleidingen.nl)

E-mail: info@sferd.nl



METHODS

Duodenoscopes

- 8 duodenoscopes with disposable endcaps

Period

- February 2022 until December 2023

Sampling



- Prior to ERCP-procedure
 - Swab distal tip
 - Flush-brush-flush sample of suction and biopsy channel (20 ml sterile water / single use endoscope pull through cleaning brush)
- Contamination = presence of microorganisms of gastrointestinal or oral origin (MGO) (including *Pseudomonas aeruginosa* and *Staphylococcus aureus*)

Data collection


- Data collection from electronic patient records on use and reprocessing of duodenoscopes and reprocessing staff



POTENTIAL RISK FACTORS

1 (Bedside) pre-cleaning  >30 minutes  Manual cleaning²

2 ≤ 5 minutes manual cleaning

3  Reprocessing employee <30 reprocessing cycles over study period

Drying/storage:

4 < 90 minutes (incomplete drying cycle)³

5 >7 days²

2. Beilenhoff et al, ESGE, Endoscopy; 2018

3. Instructions for use, drying cabinet



RESULTS



307 duodenoscope cultures



58 (18.9%) positive with MGO

Baseline characteristics reprocessing cycles, n=1296 (100%)

| | |
|--|--------------|
| Duodenoscope uses, n (%) | 1150 (88.7%) |
| Delay >30 minutes before reprocessing, n (%) | 502 (47.2%) |
| Manual cleaning ≤5 minutes, n (%) | 212 (16.9%) |
| Reprocessing employee <30 reprocessing cycles, n (%) | 308 (24.4%) |
| Drying time <90 minutes, n (%) | 61 (4.8%) |
| Storage >7 days, n (%) | 13 (1.0%) |



MIXED-EFFECT MODEL

MGO

Manual cleaning time 5 min. or less

Delay >30 min. before reprocessing

Reprocessor with <30 reprocessing cycles

Drying time <90 min.

Storage >7 days

Duodenoscope use

adjusted OR (95% CI)

P-value




DOES FREQUENCY MATTER?

| | HCW \geq 30 cycles over study period | HCW $<$ 30 cycles over study period | p-value |
|--|--|-------------------------------------|-----------|
| Reprocessing cycles performed, n | 955 (100%) | 308 (100%) | |
| Reprocessing employees, n | 6 | 35 | |
| Reprocessing cycles performed per HWC, median [Q1, Q3] | 62.50 [31.00, 237.00] | 6.00 [4.00, 17.00] | 0.001 |
| Delay $>$ 30 minutes before reprocessing, n (%) | 376 (45.5%) | 126 (52.9%) | 0.05 |
| Manual cleaning \leq 5 minutes, n (%) | 130 (13.7%) | 82 (26.9%) | $<$ 0.001 |

HCW: Health Care Worker



SENSITIVITY ANALYSIS



SUMMARY AND CONCLUSIONS

- **Duodenoscope contamination**
 - 18.9% contamination rate
 - Preventable healthcare-associated infections
- **Risk factors within reprocessing**
 - Manual cleaning ≤ 5 min.
 - Indirectly: Frequency of duodenoscope reprocessing by reprocessing staff
- **Reevaluation reprocessing protocols**
 - Minimum time manual cleaning
 - Reprocessing staff: minimum frequency of reprocessing?
- **Digital monitoring and alert systems beneficial?**



THANK YOU FOR YOUR ATTENTION



E-mail: a.bulkmans@erasmusmc.nl

